

「谷浦式 CT ポジショナー」 (Ahmics-VMD HT-S1)

四肢保定具_犬種別適合サイズ

【後肢】

	5	10	15 (kg)
チワワ			
ヨークシャ・テリア			
ポメラニアン			
マルチーズ			
シー・ズー			
トイプードル			
M・ダックス			
パグ			
柴			
ビーグル			
M・シュナウザー			
シェルティ			
雑種			
★ペキニーズ			
★フレンチ・ブルドック			
コーギー			
猫			

 小  中  大

「谷浦式 CT ポジショナー」 (Ahmics-VMD HT-S1)

四肢保定具_犬種別適合サイズ

【前肢】

	5	10	15 (kg)
チワワ	[Blue bar from 3 to 4.5]		
ヨークシャ・テリア	[Blue bar from 3 to 4.5] [Green bar from 4.5 to 5.5]		
ポメラニアン	[Blue bar from 3 to 4.5] [Green bar from 4.5 to 5.5]		
マルチーズ	[Blue bar from 3 to 4.5] [Green bar from 4.5 to 6.5]		
シー・ズー	[Blue bar from 3 to 4.5] [Green bar from 4.5 to 6.5]		
トイプードル	[Blue bar from 3 to 4.5] [Green bar from 4.5 to 6.5]		
M・ダックス	[Blue bar from 4 to 6.5]		
パグ	[Green bar from 4 to 7.5]		
柴	[Green bar from 5 to 8.5]		
ビーグル	[Green bar from 6 to 11.5] [Purple bar from 11.5 to 13.5]		
M・シュナウザー	[Green bar from 5 to 8.5]		
シェルティ	[Green bar from 7 to 11.5] [Purple bar from 11.5 to 13.5]		
雑種	[Green bar from 5 to 11.5] [Purple bar from 11.5 to 13.5]		
★ペキニーズ	[Blue bar from 4 to 6.5]		
★フレンチ・ブルドック	[Green bar from 7 to 10.5]		
コーギー	[Blue bar from 5 to 13.5]		
猫	[Green bar from 3 to 5.5]		

小
 中
 大